



Standard Colors:



Red Velvet
(Installed)



Coastal Blue
(Installed in Box)

ALL NEW Sporty Blacked-Out Styling!

Great for Indoor & Outdoor Use!

Buzzaround 3-Wheel GB107



The new and improved Buzzaround LT 3-Wheel scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Upgrades include a sleek and sporty blacked-out design, a lightweight padded stadium-style seat, contoured hand grips, and a longer adjustable tiller for the ultimate driving experience. The new Buzz LT features a color display on the console, clear-coated plastic color through shrouds (red installed and blue included in the box), scuff guard on the foot deck, and additional side reflectors for added safety. With a weight capacity of 300 lbs., it also has increased foot room, a convenient charging port on the tiller, and a rugged hard plastic basket. Enjoy its indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance!

Standard Features:

- 300 lb. weight capacity
- Lightweight, padded stadium-style seat
- Ergonomic hand grips
- Non-slip deck with scuff guard
- Clear-coated, durable color through shrouds
- Battery charging port on the tiller
- Side reflectors
- Disassembles quickly and easily for transport



FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.



LIGHTWEIGHT PORTABLE SCOOTER

Model	GB107
Medicare Code	K0800
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range ¹	7 miles
Ground Clearance Under Deck	2.5"
Turning Radius	39"
Type Batteries	2-12AH
Colors	Red & Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	50 Amp
Dimensions:	
Length	41.75"
Width	19-21"
Height (ground to top of seat back)	35" - 36.5"
Ground to Top of Seat	21" - 22.5"
Ground to Top of Deck	4"
Deck to Top of Seat	15" - 17"
Number of Seat Height Positions	3
Size of Increments	.75"
Front Axle to Rear Axle	33.5"
Weight of Unit: (assembled)	85.5 lbs
Front Half	28 lbs
Rear Half	20.5 lbs
Seat w/Arms	19 lbs
Batteries (pack with 2 batteries)	17.5 lbs
Tires:	Solid, Low Profile
Front	8"
Rear	8"
Rear Anti-Tippers	2"
Standard Seat:	
Back Height - no headrest	15"
Width x Depth	17.5"x15"
Color	Black Vinyl
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Battery	1 Year



Sporty black wheels and an impressive turning radius make tight turns effortless.



Pivoting tiller folds easily.



Sleek dash with comfort grip handles and conveniently located charging port on the tiller.



Padded, contoured seat with width-adjustable, flip up arm rests for customized comfort.



Quickly disassembles into five, lightweight pieces for easy transport!

¹. Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

². Electronics warranty excludes batteries.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Find us on
Facebook

www.facebook.com/goldentechnologies